

Connect the camera to the mac and insert the dongle (USB-like stick). Boot to windows (hold alt key) and login. **Login/pw: student**
 Double-click on the ZEN2011 desktop icon

Before taking an Image

- set file names: **Tools (3) -> Options: Naming**
- set Category to **Snap** and enter **Prefix** (name/ initials, group, course)
- add date as predefined **Format IDs**
- set chosen objective in the tool bar (**center screen area**)
- make sure **Scaling automatic** is set in the status bar




Image Acquisition


- set beamsplitter to camera (out)
- open the **Acquisition** tab
- select **4ChFluorescence** (1)
- select channels (2)
- click **Live** button, (live image is displayed in center screen area)
- adjust focus, if required
- activate **Range Indicator** (view controls area), blue = under/ red = over-saturated
- **Set Exposure** time automatically, then adjust using the **Time** slider
- click **Snap** to start the acquisition
- > follow the instructions for switching reflector and/ or light source
- open **Histo** tab (center screen area) to **check for over-saturation!**

Enhance visualisation tools

- open **view controls area** and select **Dimensions** tab:
- select visualized channel(s) and change/ remove pseudo colors
- Display** tab:
- select channel(s) to be modified
- use **Min/Max** or **Best Fit** to improve visualisation

Save/ Export your files (for ImageJ analysis in ZVI format!)

-> open **Processing** tab -> select **Batch** -> set batch method to **ZVI Export**



- select output folder or Use Input Folder as Output Folder
- files and define output- (as described above)
- select files and to export files

Use File -> Save all... (select output folder) to save in CZI format
Clean up default folder when you are done: C:\Users\student\Pictures\

